


Name _____

Date _____

EQUIPMENT FOR SEPARATION

Separation equipment used in chemistry labs includes funnel, filter flask, water bath and dropping bottle.

Different separation methods



The following are questions about the equipment used.

1. Can you match the equipment to its purpose?

SORT ELEMENTS - Filtering flask, Funnel, Dropping bottle, Magnets, Sieve, Filter paper

It is used to reduce pressure.

It is often used to separate solid and liquid mixtures.

It is used to add liquid by drop number, always with a rubber head.

It is used to sift earth materials.

It is used to separate non-magnetic material from magnetic material.

It is used to sieve solids at different phrases.